

DipSensor[®]

The Farmer's Friend

Manage antibiotic risk in milk
with confidence

*Easy sampling anywhere
Your result within 10 minutes
Easy interpretation*

For
farm
use

Have you ever had concerns that your milk might be contaminated?

Has a cow lost its tag, has the cow calved early or is a cow still antibiotic positive after the withdrawal period?



DipSensor helps you to manage the antibiotic risk at your farm

If the antimicrobial risk within a herd is not evaluated properly, the consequences for a farmer can be extremely serious:

- Milk tests positive for antibiotics by the dairy, and resulting in penalties
- Contamination of the whole milk collection truck

Do you want to be sure? DipSensor gives you that confidence right at the farm!

- No equipment needed
- Easy sampling anywhere
- Your result within 10 minutes
- Easy interpretation


DipSensor®



Place the Dipper
into the milk




Insert the test strip
in the Dipper
in contact with the milk



Incubate
for 10 minutes
at room temperature

Visual Interpretation




Sensitivity table

BETA-LACTAMS						TETRACYCLINES		
PENICILLINS			CEPHALOSPORINS					
	EU MRL (PPB)	LOD (PPB)		EU MRL (PPB)	LOD (PPB)		EU MRL (PPB)	LOD (PPB)
Amoxicillin	4	3-4	Cefacetrile	125	15-30	4-Epimer of chlortetracycline	100	60-100
Ampicillin	4	2-4	Cefalexin	100	50-100	4-Epimer of oxytetracycline	100	60-100
Aspoxicillin	-	5-10	Cefalonium	20	2-5	4-Epimer of tetracycline	100	60-100
Cloxacillin	30	10-15	Cefazolin	50	10-20	Chlortetracycline	100	60-100
Dicloxacillin	30	5-10	Cefoperazone	50	1-4	Doxycycline	-	30-60
Nafcillin	30	80-150	Cefquinome	20	20-40	Oxytetracycline	100	30-60
Oxacillin	30	15-20	Ceftiofur	100*	20-40	Tetracycline	100	60-100
Penicillin G / Benzylpenicillin	4	2-3	Desfuroyl ceftiofur		200-300	* The MRL is applied to the sum of both residues		
Phenoxymethylpenicillin	-	2-4	Cefuroxime	-	200-300			
Piperacillin	-	1-4	Cephapirin	60*	4-8			
Ticarcillin	-	15-30	Desacetylcephapirin		30-50			

KIT900x25

Each bag includes:

- 25 test strips
- 5 dipper tubes



No instrument needed!
Room temperature
10 minutes



Unisensor - Who We Are

Our Mission

We bring smart diagnostic solutions from the laboratory to the field.

Throughout our history, we have responded to the needs of the field with new technologies that are unmatched in the market.

Through innovation, we don't just optimise our products within our business, we aim to reshape and redefine our industry.

Our Vision: Peace Of Mind in Food

Consumer confidence is essential for food producers. With our products, we help them guarantee superior quality food that meets consumers' expectations.





Liège Science Park • Rue Louis Plescia, 8 • 4102 Ougrée (Liège) • Belgium
Phone +32 4 252 66 02 • info@unisensor.be • www.unisensor.be



Chaussée de Wavre 1110 - 1160 Bruxelles
Waversesteenweg 1110 - 1160 Brussel
02/660.50.75 - info@intermed.be - www.intermed.be