# IC BOVINO IC CAPRINO

ZE/ICB

ZE/ICB100

ZE/ICC

ZE/ICC100

**Immunochromatographic test for  
milk species identification**

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info@zeulab.com • www.zeulab.com

**SCOPE**

IC BOVINO kit is a qualitative test to detect the presence of cow's milk in goat or sheep's milk. The test is an immunoassay based on the detection of bovine immunoglobulins (IgG) in milk and cheese samples.

IC CAPRINO kit is a qualitative test to detect the presence of goat's milk in sheep's milk. The test is an immunoassay based on the detection of caprine immunoglobulins (IgG) in milk and cheese samples.

<b>KIT COMPONENTS</b>	<b>ZE/ICB &amp; ZE/ICC</b>	<b>ZE/ICB100 &amp; ZE/ICC100</b>
<i>Strip tests</i>	25	100
<i>Dilution solution</i>	1	2
<i>Disposable pipettes</i>	25	100
<i>Test tubes</i>	25	-
<i>Divisible strip of 8 wells</i>	-	13
<i>Tube/well rack</i>	yes	yes

**SAMPLES PREPARATION AND SAMPLE INFORMATION**

Whole, skimmed, fresh and pasteurised milk; whey samples and cheese can be used with these kits. Nevertheless, be aware that heat treatments over 72°C/15-20 sec can begin to denaturize milk immunoglobulins. The higher the heat treatment, the higher the immunoglobulins denaturalization will be. Therefore, if possible, test samples prior to being processed.

Milk samples can be stored under refrigeration for 1 or 2 days before the assay. For longer period of time, samples should be frozen at -20°C.

**Milk sample preparation:** No sample preparation is needed

**Cheese sample preparation:**

1. Weigh 5g of an homogeneous cheese sample and add 10 mL of distilled water and mix with a pestle and mortar to obtain a fine homogenised sample.
2. Centrifuge at 3000 x g for 10 minutes. Alternatively, the sample can be filtered using a filter paper.
3. The supernatant or the filtered solution will be used for testing.

## ASSAY PROCEDURE

1. Take a test tube or a well from the divisible strip. Place it into the rack provided.
2. Dilute the sample using a new disposable pipette for each sample. Stir gently.
  - **Milk samples:** mix 3 drops of dilution solution with 1 drop of milk sample.
  - **Cheese samples:** mix 2 drops of dilution solution with 1 drop of the cheese extract.
3. Place the strip test into the test tube/well following the arrows direction.
4. Wait for 10 minutes and read the results within the next 2 minutes.
5. Read the result in the following 2 minutes after the end of the analysis.

### Visual reading

**NEGATIVE:** A BLUE line is observed in the white central zone of the strip test. This line indicates that the test has been run properly.

**POSITIVE:** Two lines, BLUE and RED, are observed in the white central zone of the strip test.

### Objective reading - IC-Reader

Place the strip in the adapter and couple it with the IC-Reader. Place the RFID card corresponding to the kit batch on top of the device and press the button. Please, check IC-Reader instructions for further information.

Results will be displayed on the screen as follows:

-Milk samples: **NEG** (negative), **POS** (positive) or **HPOS** (high positive).

For more information about levels of mixtures, please check the product certificate.

-For cheese: **NEG** (negative) or **POS** (positive).

Note: There is a specific reading method for each batch and for the analysis of milk or cheese. The method is transferred to the reader memory via RFID card and is still available after the IC-Reader is turned off. When different batches are used, pay special attention to the method uploaded on the reader.

**NOTES**

- The kit is designed to detect traces of milk mixture. As example, 0.125% of mixture implies that a tank of 10 000 L is contaminated with 12.5 L of other milk specie.
- The red line of the test might not appear when testing samples of pure milk or mixtures of more than 70% of the target milk due to "prozone or Hook effect".
- The intensity of the line might change over time. Results should be read in the following 2 minutes after the end of the analysis.
- To keep records of the result use IC-Reader or take pictures of the strips when the assay is finished.
- If the blue line does not appear on the white central zone, the test has not been run properly or the reagents have been damaged. The result must be considered INVALID (INV) and repeat it. If it happens again, please contact ZEULAB.

**DETECTION LIMITS**

Limit of Detection in milk: IC Bovino 0.25% of mixture. IC Caprino 0.5% of mixture.

**STABILITY**

Store the test under dry conditions and temperatures between 15 and 25°C. The strips are very sensitive to humidity and must be kept inside the container provided with the desiccant bag. The kit has a shelf life of approximately 12 months. Check the expiration date on the package.

**SAFETY**

This kit should be used following good laboratory practices. A Material Safety Data Sheet is available on request.

IC KIT is an *in vitro* diagnostic test. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility.

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