

Rapid Test

Detecting Tetracycline

In Any Animal Tissues




Tetracyclines are broad-spectrum antibiotics that are extensively used in veterinary and human medicine for the treatment of a wide variety of bacterial infections.


TetraSensor is a rapid test to detect the contamination of meat samples by Tetracycline molecules.

The test can be used in-situ for daily controls or in laboratories for the analysis of series of tissue samples.


BENEFITS

 **NO NEED TO HEAT / ROOM TEMPERATURE**

 **RAPID**
10 minutes incubation only.

 **SENSITIVE**
Low limits of detection.


 **RELIABLE AND ROBUST**

 **GENERIC**
Detection of all relevant Tetracycline molecules in one single test.

 **USER-FRIENDLY**
Easily performed on site or in the lab.

RESULTS

Visual interpretation of the result obtained by comparing the intensity of each test line with the control (CTRL) line.



An instrumental reading is also possible with the ReadSensor 2 device.

SENSITIVITY


| COMPOUNDS | LIMITS OF DETECTION (PPB) | |
|-------------------|---------------------------|-------------------|
| | CRUDE TISSUE (4X)* | COOKED MEAT (6X)* |
| Chlortetracycline | 2 - 3 | ND |
| Doxycycline | 2 - 3 | 5 - 10 |
| Oxytetracycline | 5 | ND |
| Tetracycline | 10 | ND |

* MATRICES (Considering dilution)

HOW TO USE

SAMPLE PREPARATION


ROOM
TEMPERATURE




Full procedure available for sample preparation of :

- Raw muscle tissue
- Processed tissue (including cooked meat)

TEST EXECUTION



Place the vial to plastic tray




Mix sample
and reagents



Add dipstick
to the vial

INCUBATE **10** MINUTES
AT ROOM TEMPERATURE



Interpret
results

PRODUCTS REFERENCES

KIT036
TETRASENSOR TISSUE

APP088
READSENSOR 2



unisensor
DIAGNOSTIC ENGINEERING

Liège Science Park
Rue Louis Plescia, 8
B-4102 Ougrée (Seraing)
BELGIUM

www.unisensor.be

Phone +32 4 252 66 02
info@unisensor.be

Reseller

im NV International Medical Products SA
Marcel Broodthoersplein 8/5 Place Marcel Broodthoers
BRUSSEL 1060 BRUXELLES - BELGIUM
T: +32 (0)2 660 50 75 - F: +32 (0)2 660 20 98
www.intermed.be - info@intermed.be

beldico B.V.
Verlengde Poolseweg 16
4818 CL Breda
T: +31 (0) 765 244 663
www.beldico.nl - info@beldico.nl