

4Sensor BSCT



The test can be used in-situ for daily controls or in laboratories for the analysis of series of milk samples. This test does not require any sample processing, cleaning or extraction.

BENEFITS



MULTIPLE

Detection of Betalactams, Tetracyclines, (Dihydro)streptomycin and Chloramphenicol.



RAPID

Results in 10 minutes.



SENSITIVE

EU MRLS.



RELIABLE AND ROBUST



USER-FRIENDLY


Easily performed on site or in lab.



COST-EFFECTIVE

RESULTS

Visual interpretation of the result obtained by comparing the intensity of each test line with the control (CTRL) line.



T: Tetracyclines - C: Chloramphenicol - S: Streptomycin - B: Betalactams

An instrumental reading is also possible with the ReadSensor 2 device.

MEETS REQUIREMENTS OF THE CUSTOMS UNION





The development of the test was made in connection with Danone-Russia and All-Russia Dairy Industry Research Institute.


GOST R53774-2010 APPROVED GOST 32219-2013



SENSITIVITY

	ß-LACTAMS	LOD (PPB)
PENICILLINS	Penicillin-G	2 - 3
	Ampicillin	3 - 4
	Amoxicillin	
	Oxacillin	12 - 18
	Cloxacillin	
	Dicloxacillin	6 - 8
	Nafcillin	30 - 40
	Ceftiofur	10 - 15
	Cefquinome	30 - 35
	Cefazolin	18 - 22
	Cefapirin	
	Cefacetrile	30 - 40
	Cefoperazone	
	Cefalexin	1000 - 1200
	Cefalonium	3 - 5
	TETRACYCLINES	LOD (PPB)
	Tetracycline	8 - 10
	Oxytetracycline	7 - 9
	Chlortetracycline	
	Doxycycline	5 - 7 2 - 3
	Doxycycline	2 - 3
	Doxycycline CHLORAMPHENICOL	2 - 3 LOD (PPB)
	Doxycycline CHLORAMPHENICOL Chloramphenicol	2 - 3 LOD (PPB) 0,3

HOW TO USE



Mix milk and reagents



Place the vial to heating block




Close and push "START"



Place dipstick to Heatsensor





4Sensor BSCT

PRODUCTS REFERENCES

KIT060

4 SENSOR MK MRL - 96 PC

APP032
HEATSENSOR DUO
APP052
HEATSENSOR OCTO



APP088 READSENSOR 2





Liège Science Park Rue Louis Plescia, 8 B-4102 Ougrée (Seraing) BELGIUM

www.unisensor.be

Phone +32 4 252 66 02 info@unisensor.be

