

## Quantitative Rapid Test Detecting

## Aflatoxin M1 In Milk

AflaSensor is a rapid assay in dipstick format detecting quantitatively the contamination of milk samples by Aflatoxin M1 molecules.

The test can be used in-situ for daily controls or in laboratories for the analysis of series of milk samples. This test does not require any sample processing, cleaning or extraction.



## BENEFITS

**QUANTITATIVE****RAPID**

Results in 10 minutes.

**SENSITIVE**

Adapted to specific regulation (EU and US).

**RELIABLE AND ROBUST****USER-FRIENDLY**

Easily performed on site or in lab.

**COST-EFFECTIVE**

## RESULTS

The use of Readsensor 2 allows an accurate quantification of the Aflatoxins contamination.



## SENSITIVITY


TOXIN	EU MRL	US ACTION LEVEL (SL)
	KIT041 (QUANTITATIVE)	KIT078 (QUANTITATIVE)

Aflatoxin M1 (PPT)


Detection range from 30 to 150

Detection range from 200 to 750

## HOW TO USE




1  
Mix milk  
and reagents



2  
Place the vial  
to heating block



3  
Close and  
push "START"



4  
Place dipstick  
to Heatsensor

INCUBATE **10** MINUTES



5  
Interpret  
results

QUANTITATIVE

# AflaSensor Quanti

## PRODUCTS REFERENCES

### KIT041v2

AFLASENSOR QUANTI MK 0,05 PPB - 96 PC

### KIT078

AFLASENSOR MK 0,5 PPB 96 pc

APP032 HEATSENSOR DUO

APP008 AIRSHIELD DUO

APP052 HEATSENSOR OCTO

APP033 AIRSHIELD OCTO



APP088

READSENSOR 2



Liège Science Park  
Rue Louis Pléscia, 8  
B-4102 Ougrée (Liège)  
BELGIUM

[www.unisensor.be](http://www.unisensor.be)

Phone +32 4 335 30 00  
[info@unisensor.be](mailto:info@unisensor.be)

## Reseller

 **International Medical Products SA**  
Marcel Broodthaersplein 8/5 Place Marcel Broodthaers  
BRUSSEL 1060 BRUXELLES - BELGIUM  
T: +32 (0)2 660 50 75 - F: +32 (0)2 660 20 98  
[www.intermed.be](http://www.intermed.be) - [info@intermed.be](mailto:info@intermed.be)

 **beldico B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda  
T: +31 (0) 765 244 663  
[www.beldico.nl](http://www.beldico.nl) - [info@beldico.nl](mailto:info@beldico.nl)