

HeatSensor

HEATSENSOR® DUO - OCTO are incubators ensuring all Unisensor diagnostic tests to receive a fixed incubation temperature no matter what the ambient / operating temperature.

Both incubators are a **One Step Test Procedure**.

KEY-POINTS

- „ Automatic release of dipstick.
- „ Compatibility with all test kits from UNISENSOR.
- „ DUO is for maximum 2 tests at a time.
- „ OCTO is for maximum 8 tests at a time.
- „ Self-timing with light and audio signal.
- „ Countdown timer display.
- „ Product's specific program included.
- „ Easy to use.



SPECIFIC PROGRAM

Select the right program and place the dipstick.



AUTOMATIC

The dipstick falls down automatically into the microwell after the first incubation time.

PRODUCTS REFERENCES

APP032
HEATSENSOR DUO

APP052
HEATSENSOR OCTO



Unisensor s.a.

Liège Science Park
Rue Louis Plescia, 8
B-4102 Ougrée (Seraing)
BELGIUM

www.unisensor.be

Phone +32 4 252 66 02
info@unisensor.be

Reseller

im NV International Medical Products s.a.
Marcel Broodthaersplein 8/5 Place Marcel Broodthaers
BRUSSEL 1060 BRUXELLES - BELGIUM
T: +32 (0)2 660 50 75 - F: +32 (0)2 660 20 98
www.intermed.be - info@intermed.be

beldico b.v.
Verlengde Poolseweg 16
4818 CL Breda
T: +31 (0) 765 244 663
www.beldico.nl - info@beldico.nl