

# HeatSensor

HEATSENSOR® DUO - OCTO are incubators ensuring all Unisensor diagnostic tests to receive a fixed incubation temperature no matter what the ambient / operating temperature.

Both incubators are a **One Step Test Procedure**.

## KEY-POINTS

- Automatic release of dipstick.
- Compatibility with all test kits from UNISENSOR.
- DUO is for maximum 2 tests at a time.
- OCTO is for maximum 8 tests at a time.
- Self-timing with light and audio signal.
- Countdown timer display.
- Product's specific program included.
- Easy to use.



### SPECIFIC PROGRAM

Select the right program and place the dipstick.



### AUTOMATIC

The dipstick falls down automatically into the microwell after the first incubation time.

## PRODUCTS REFERENCES

**APP032**  
HEATSENSOR DUO

**APP052**  
HEATSENSOR OCTO



Unisensor S.A.

Liège Science Park  
Allée de la Cense-Rouge, 98  
B-4102 Ougrée (Liège)  
BELGIUM

[www.unisensor.be](http://www.unisensor.be)

Phone +32 4 252 66 02  
Fax +32 4 252 90 55  
[info@unisensor.be](mailto:info@unisensor.be)

### Reseller

 **N.V. International Medical Products S.A.**  
Waverssesteenweg 1110/7 Chaussée de Wavre  
OUDERGEM 1160 AUDERGHEM - BELGIUM  
T: +32 (0)2 663 15 45 - F: +32 (0)2 660 20 98  
[www.intermed.be](http://www.intermed.be) - [info@intermed.be](mailto:info@intermed.be)