

# Maximize AFB Recovery, Minimize pH-Related False Positives

## A Unique Digestion and Decontamination System

*for the recovery of Mycobacterium spp.*

### The Alpha-Tec NAC-PAC™ RED Advantage

The essential first component in the Alpha-Tec mycobacteria specimen preparation system

Provides an effective, visual confirmation of a neutral pH without titration

The NALC and NAC-PAC RED mixture supports mucolytic activity for up to 72 hours

Available in multiple NaOH concentrations and volume offerings

### The NAC-PAC RED System

Each NAC-PAC RED package contains five pre-measured NALC ampules and five bottles of pre-mixed and filtered TB base digestion solution with pH indicator, unless otherwise noted

Product Code	Product Name	Quantity
0004302	NAC-PAC RED, 2.0%	5 x 50 ml
0004307		5 x 200 ml
0004305	NAC-PAC RED, 2.5%	5 x 50 ml
0004306		5 x 200 ml
0004310	NAC-PAC RED, 3.0%	1 x 200 ml No NALC ampules
0004303		5 x 50 ml
0004309		5 x 200 ml
0004311		1 x 200 ml No NALC ampules
0004304	NAC-PAC RED, 4.0%	5 x 50 ml
0004308		5 x 200 ml

# NACPAC™ RED

