



COVID-19 Ag K-Set

Detection of SARS-CoV-2 antigen
in nasal secretions

HIGH SENSITIVITY, EASY AND RAPID




Performances


Compared to RT-PCR (n=272)

		95 % Confidence Interval ¹
Sensitivity	96.9 %	(90.7 to 99.2 %)
Specificity	100 %	(97.3 to 100 %)
Positive predictive value	100 %	(95.2 to 100 %)
Negative predictive value	98.3 %	(94.7 to 99.6 %)
Agreement	98.9 %	(269/272)

Procedure



Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

 **International Medical Products SA**
Waverssesteenweg 1110/7 Chaussée de Wavre
OUDERGEM 1160 AUDERGHEM - BELGIUM
T: +32 (0)2 660 50 75 - F: +32 (0)2 660 20 98
www.intermed.be - info@intermed.be