


GastroVir K-Set


Rapid diagnostic test for *in vitro* detection of Rotavirus & Adenovirus (serotypes 40/41) in stool specimens

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (10')
- Multiple coloured lines



www.corisbio.com

GastroVir K-Set



Performance of the kit


	Rotavirus detection Compared to PAGE method (n=400) (95% confidence interval)	Adenovirus 40/41 detection Compared to PCR method (n=400) (95% confidence interval)
Sensitivity	97.4% (92.1 to 99.3%)	96.7% (80.1 to 99.8%)
Specificity	97.9% (95.2 to 99.1%)	97.3% (94.9 to 98.6%)
PPV	95.0% (88.9 to 97.9%)	74.4% (57.6 to 86.4%)
NPV	98.9% (96.7 to 99.7%)	99.7% (98.2 to 100%)
Agreement	97.8% (391/400)	97.3% (389/400)

For more technical information

client.care@corisbio.com

Ordering Information

PRODUCT NAME	DESCRIPTION	STORAGE	CODE
GastroVir K-Set with collection set	20 tests, 20 faecal sampling systems	4°C to 30°C	K-1516



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

Distributed by:



D-CUS-10-16-K/07
March 2022