




VitaPCR™

Easy. Rapid. Versatile.





POINT-OF-CARE MOLECULAR ASSAY

FEARURES & BENEFITS

✔ Easy

- Intuitive: only 2min hands-on time
- Minimal training required
- Automated result interpretation

✔ Rapid

- Sample to result in less than 20min
- The right diagnosis on the spot
- Improved patient management

✔ Versatile

- Can be deployed anywhere
- No additional equipment required
- Ready-to-use reagents




2 minutes hands-on time

SPECS SARS-CoV-2


Sample Collection	NSP or OP swab
Sensitivity	2,73 copies/ μ l
Dimensions	205 x 165 x 155mm
Weight	1,2kg

EASY TO USE AND FAST




1. COLLECT

NSP or OP swab sample collection



2. LYSE

Stir the swab in the collection buffer



3. TRANSFER

Use the disposable pipette




4. START

Results delivered in less than 20 minutes


TOTAL PROCEDURE IN LESS THAN 20 MINUTES

IMPACT ON PATIENT MANAGEMENT

Laboratory Batch PCR



VitaPCR





PROVIDING ADVANCED MEDTECH FOR ALL.



BELDICO BV

Verlengde Poolseweg 16 - 4818 CL Breda

0765 244 663 - info@beldcio.nl - www.beldico.nl