

# **iPeak**®

# RAPID SEROLOGICAL TEST READER QUICKZEN COVID-19 IgM/IgG



## iPeak® Reader

The **iPeak® Reader** is a **rapid test cassette reader** that provides accurate and reliable results with high traceability. This reader can objectivize the reading of quick tests for the detection of **COVID-19** (**SARS-COV-2**) antibodies directly at the end of the ZenTech kit's incubation time: **QuickZen Covid-19 IgM/IgG**.

The iPeak® Reader is compact and user-friendly. The embedded software enables the objective and qualitative interpretation of the test. The results can then be exported and data integrity and traceability are guaranteed.



Reproducibility and precision



Integrated barcode reader (ID – 2D – QR Code)



**Data traceability** 



# **Device description**

The iPeak® Reader is a portable easy-to-use system. It provides quick and accurate results with a high level of sensitivity. When used with the 'QuickZen COVID-19 IgM/IgG' ZenTech kit, the iPeak® Reader can read test lines and applies a special algorithm to interpret results for each of these lines.

With its 7 inch capacitive touchscreen, the iPeak® Reader is the ideal tool for taking and interpreting test images on samples in real time.

It features QR Code/Datamatrix reading for full sample traceability and a range of connection options for exporting the results.



A **printer** can also be connected to the **USB port**.

The iPeak® Reader can store a large volume of results, raw images and data in its internal memory. It works on **battery** and/or **mains**.

Reference: H-LY-000

#### **Performance**



Results displayed within a minute



**Internal battery** (up to 8h autonomy)



Results stored in the system's internal memory



7 inch touchscreen



Data exported via Ethernet, Wi-Fi and USB port Wifi



**Printer connection option** 

### Related kit



Packaging: 25 individual tests

Reference: S-LW-25

Kit content: 25 pouches containing one cassette, 25 whole blood collection pipettes, a 4ml vial of dilution buffer, one instruction

manual.



Rapid 10 minutes



**Practical** 



Reliable All-in-one Excellent sensitivity and specificity



#### ZenTech S.A.

Liège Science Park Avenue du Pré-Aily, 10 - 4031 ANGLEUR - BELGIUM Phone: +32 (0)4 361 42 32 - Fax: +32 (0)4 367 00 63

info@zentech.be - www.zentech.be

#### **Distributor for Belgium**

#### **International Medical Products**

Chaussée de Wavre 1110 / 7 - 1160 Brussels - Belgium Sales and customer service: + 32 (0)2 661 33 00 Contact: corona@intermed.be