Discover the whole RESIST range



| Product | Target | Product | Target |
|--|---------------------------|---------------------------------------|---|
| OXA-48 K-SeT Code: K-15R1 | OXA-48-like | OXA-23 K-SeT Code: K-15R7 | OXA-23-like |
| KPC K-SeT Code: K-15R2 | KPC | IMP K-SeT Code: K-15R10 | IMP |
| RESIST-3 O.K.N. <i>K</i> -SeT Code: K-15R5 | OXA-48-like KPC NDM | RESIST-3 O.O.K. K-SeT Code: K-15R4 | OXA-48-like OXA-163 KPC |
| | | RESIST-5 O.O.K.N.V. Code: K-15R9 | OXA-48-like OXA-163 KPC NDM VIM |

- Content of the kit: cassettes individually pouched, buffer, tubes with droppers, Instruction for Use
- **Storage:** 4°C to 30°C


CORIS

- Stability: from 15 to 24 months at RT
- Bibliographic and Technical support: client.care@corisbio.com
- Price: please refer to your local distributor



Solution adapted to your daily and epidemic needs for the detection and identification of carbapenemases.

O.K.N.V.I. RESIST-5








Science Park CREALYS Rue Jean Sonet 4A - 5032 Gembloux - Belgium Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919 sales@corisbio.com - www.corisbio.com



D-CUS-10R11/02

OXA-48 - KPC - NDM - VIM - IMP









Choose the right tool for your carbapenemase detections. We offer you the kit adapted to your day-to-day needs.

RESIST OXA-48

RESIST OXA-23

RESIST KPC

RESIST IMP

RESIST-3 O.K.N. = OXA-48, KPC, NDM

RESIST-3 O.O.K. = OXA-48, OXA-1 63, KPC

RESIST-5 O.O.K.N.V. = OXA-48, OXA-1 63, KPC, NDM, VIM

O.K.N.V.I. RESIST-5 = OXA-48, KPC, NDM, VIM, IMP

