

Bottle warmer & Pasteurizer Inserts



To make it easier for nursing staff to prepare their little patients' meals, Beldico has developed inserts so that different sizes of feeding bottles and enteral syringes can be warmed in the same cycle.

Fitting perfectly into the cavities of the Beldico bottle warmers and pasteurizer, and made from the same material as these devices' heating blocks, the inserts ensure that heat is evenly and constantly distributed throughout the milk.

The inserts are replaceable and can be cleaned and disinfected by hand or washed in a dishwasher*. Their smooth surface makes them easy to wash or sterilise, to ensure hygienic warming or pasteurization of the milk prepared for the baby.

* Contact us for a list of disinfection products that are compatible with Beldico bottle warmers and inserts.

Adapters for small Beldico bottle sizes

8618007	AE-RED-4	50 ml bottle adapter for 1-2-4-8-BA120 / 2-4-8-BA240**
5101891	PA45-50ML	50 ml bottle adapter for PA45
5100112	AE-RED-6	130 ml bottle adapter for 2-4-8-BA240 / PA45
5100113	AE-RED-6 stick	Stick for removing AE-RED-6 adapters

** Use AE-RED-4 in combination with AE-RED-6 for 2-4-8-BA240 bottle warmers

Inserts for other special bottles

8618008	AE-RED-1	Insert for Nestlé Pre NAN 90 ml
8618005	AE-RED-2	Normal insert for Friso 50 ml
8618006	AE-RED-3	Standard insert for Nutrilon 60 ml Insert for Milupa Aptamil Primergen 60 ml Insert for Nutrilon Neonatal Start 60 ml
5100159	AE-RED-7	Insert for Milupa Aptamil Malto 15 100 ml
5100160	AE-RED-8	Insert for Humana Anfangsmilch PRE 90 ml Insert for Nutrilon Hypo-Allergen 90 ml
5100397	AE-RED-9	Insert for Milupa Nutrilon Bledina (Plastic) 90 ml



List not exhaustive, consult us for other formats.

Inserts for enteral syringes

Insert for 1 ml syringes – 5 cavities

Insert for 2 ml syringes – 4 cavities

Insert for 5 ml syringes – 3 cavities

Insert for 10 ml syringes – 3 cavities

Insert for 20 ml syringe – 1 cavity

Insert for 35 ml syringe – 1 cavity

Insert for 60 ml syringe – 1 cavity


Insert for 2.5 ml + 5 ml+ 10 ml syringes – 3 cavities

Insert for 2 ml + 2 x 10 ml syringes – 3 cavities

Insert for 5 ml + 10 ml syringes – 2 cavities

List not exhaustive, consult us for other formats, giving details of:

- Brand of syringes
- Connection type (ENFit, ...)
- Connection position (central, off-centre)
- Size of syringes used
- Use with or without an outer label on the syringe



Beldico develops, produces and markets medical devices for feeding babies and focuses particularly on the specific needs of newborn and prematurely born babies

Bottles – breastshields – teats – pacifiers
Breast milk expression and handling
Breast milk or formula warming

Manufacturer

AERNE ANALYTIC E.K.

GERMANY

www.aerne-analytic.de

Distributors

INTERNATIONAL MEDICAL PRODUCTS S.A.

Chaussée de Wavre 1110/7

BE-1160 Auderghem

Bruxelles

BELGIUM

Tel : 0032 (0) 2 660 50 75

www.intermed.be

BELDICO B.V.

Marketing 22b

NL-6921 RE Duiven

NEDERLAND

Tel : 0031 (0) 263 210 210

www.beldico.nl

BELDICO france S.A.S

ZA Node Park Touraine

Rue Gilles de Gennes-BP8

F-37310 Tauxigny

FRANCE

Tel : 0033 (0) 247 255 650

www.beldico.fr