




Duckworth & Kent Ltd.
titanium surgical instrument manufacturer

at the leading edge

NEWSLETTER

FEBRUARY 2018



It's never been easier to connect and interact with Duckworth & Kent.
Stay up to date with all our latest innovations in ophthalmology via social media and our website.

www.duckworth-and-kent.com

Find us on




CORE VALUES *Collection*

Core Values Collection

Innovation has been at the heart of our philosophy for over half a century, driven by an uncompromising quest for perfection and an enduring commitment to exceptional engineering. Duckworth & Kent is proud to present the "Core Values Collection".

6-496-2

Drysdale Rotator



- 0.8mm x 2.4mm large paddle
- 45° angled shaft, tip to angle length 10mm
- Round handle, length 118mm

Designed to be used to rotate and manipulate the nucleus.


6-656 Range

DK Lacrimal Probes

- High quality surface finish over first 5mm of tips
- Overall length 130mm




6-656




Tip sizes - 0000 & 000

6-656-1




Tip sizes - 00 & 0

6-656-2




Tip sizes - 1 & 2

6-656-3




Tip sizes - 3 & 4

6-656-4



Tip sizes - 5 & 6

6-656-5



Tip sizes - 7 & 8

For more information visit
www.duckworth-and-kent.com
to view our complete collection
of Core Values instruments.

Alternatively scan
the QR code to view
our Core leaflet.



www.duckworth-and-kent.com

Find detailed information for over 800 instruments. Our high quality reusable titanium instruments cover a multitude of ophthalmic procedures.


View accompanying literature, including our suggested sets which cover over 30 procedures: Essential Cataract, LASIK, Glaucoma, LRI, Cross-linking, Consulting Room, Chalazion, Paediatric Cataract, Vitreoretinal, Lids, Lacrimal and many more.

Stay up to date with our latest innovations in ophthalmic instrumentation.




Scleral Depressors

6-615 Tsukahara Scleral Depressor




- Highly polished tip
- Tip length 6mm, maximum diameter 1.2mm
- 55° angled shaft
- Large diameter round handle, length 114mm



Used for smoothing sclera after removing an instrument cannula used during Vitreous surgery.

6-615-1 Tsukahara Scleral Depressor with Slotted Sides




- Textured tip
- 55° angled shaft
- 3 slots on the right and left of the tip (6 slots overall)
- Tip length 6mm, maximum diameter 1.2mm
- Large diameter round handle, length 114mm



Cannulas


We have an extensive range of cannulas. See below for a selection and our website for more.

8-640 Step Instrument Cannula



- Cannula for 25 gauge Instruments
- 23 gauge version available - 8-640-2
- 45° bevelled tip, overall length 6mm
- Step prevents accidental removal

8-641-9 Infusion Cannula




- 23 gauge tip
- 45° bevelled tip
- Cannula inserted through instrument cannula
- Supplied as cannula tip only
- Overall length 6.3mm

8-641 Infusion Cannula



- 25 gauge tip only
- Cannula inserted through instrument cannula
- Silicone tubing and luer fitting supplied

8-642 Cannula Plug



- Seals 25 gauge Instrument Cannula
- Seals 23 gauge Instrument Cannula - 8-642-1
- Overall length 6.6mm
- Overall length 6.6mm - 8-642-1

View our Vitreoretinal leaflet to view a larger selection of cannulas.



E-RANGE

Instruments

E-Range

Innovation is at the heart of our philosophy and is driven by our enduring commitment to provide surgeons with only the highest quality reusable titanium ophthalmic instrumentation.

An Innovation in Ergonomic Design



This range represents the application of cutting edge technology fused with over 50 years of engineering knowledge. These instruments offer surgeons an elegant, well balanced and ergonomically designed handle whilst still providing our high quality, consistent and durable functionality.

Innovative

Ergonomic

Elegant



Capsulorhexis

Short Flat
Short Round
Long Flat
Long Round

Long Handle

Notched
Suturing
Tying
Utility

Short Handle

Notched
Suturing
Tying
Conjunctival
Cilia

Visit www.duckworth-and-kent.com for detailed information on all our E-Range instruments.

If you are interested in any of our product range or have a general enquiry please contact anyone of our dedicated customer service team members who are ready to help.



Maria
Customer Services
Supervisor



Kelly
Customer Services



Nicky

mariagente@
kellyobrien@
nickyguthrie@

duckworth-and-kent.co.uk

Find us on



© February 2018 Duckworth & Kent

Email: info@duckworth-and-kent.co.uk

Tel: +44 (0)1462 893254

Fax: +44 (0)1462 896288

www.duckworth-and-kent.com

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.