

Innovative solutions for effective diagnostics



NEW PRODUCT



Legionella K-Set

Rapid diagnostic test for *in vitro* detection of
Legionella pneumophila antigen in urine.

- One Step !
No dilution required !
- High sensitivity : 97.6%
- High specificity : 100%
- One of the best *Legionella* kit on the market !



Legionella K-SeT



● Performances ●


EIA (discrepancies classified according to EWGLI criteria n=109)		
Specificity	100%	(93.3 to 100%)
Sensitivity	97.6%	(85.6 to 99.9%)
PPV	100%	(89.1 to 100%)
NPV	98.5%	(91.1 to 99.9%)
Reliability	99.1%	

The evaluation of the new **Legionella K-SeT** kit was conducted in an European Reference Laboratory with specimens typed according to the European Working Group for Legionella Infections (EWGLI) criteria. The **Legionella K-SeT** is a one-step test whose performances rank among the best worldwide.

● Ordering Information ●

Product Name	Description	Storage	Code
<i>Legionella K-SeT</i>	20 tests, 20 transfer pipettes	4°C to 30°C	K-1215
<i>Legionella K-SeT</i>	20 tests, 20 transfer pipettes, positive control, negative control	4°C to 30°C	K-1515

Science Park CREALYS
Rue Jean Sonet 4A
5032 Gembloux
BELGIUM



Ph.: +32 (0)81 719 910
Fx.: +32 (0)81 719 919
sales@corisbio.com
www.corisbio.com