

Medimail CANISTER MCM1



PRODUCT

Sistema de TRIPLE EMBALAJE para SUSTANCIAS BIOLÓGICAS CATEGORÍA B (asignadas al número UN 3373).

//REGULATIONS//////

Packaging instruction P650 of the ADR (European Agreement Concerning International Transport of Dangerous

International Transport of Dangerous Goods by Road) and Packaging Instruction 650 de IATA/ICAO.

Triple Packaging System comprising:



Sturdy Outer Packaging made of fully-recyclable polypropylene co-polymer; resistant to water, ripping, impacts and chemical agents. It measures $290 \times 195 \times 105$ mm and can hold up to three secondary containers.

It is yellow for easy identification, with all required transport pictograms printed on it; closure is via self-folding flaps and it includes a tamperproof safety label which quarantees the security of the shipment.



A Rigid Secondary Container: this is a cylindrical container made of plastic, with a screw top and leakproof gasket to guarantee it can withstand 95 kPa of internal pressure without leakage. Its dimensions are 160 x 80 mm,(135 x 59 mm usable space),

enabling it to hold test tubes



An absorbent safety pouch made of technological polymer able to absorb up to 25 times its weight, in order to hold the total amount of liquid transported, whether bodily fluids or chemical products; it measures 170 x 120 mm and can hold up to 5 test tubes.



36208 Vigo - Pontevedra (Spain)



Medimail CANISTER MCM1

Certification.

- Complies with Packaging Instruction P650 of the ADR (European Agreement Concerning International Transport of Dangerous Goods by Road).
- Also complies with Packaging Instruction 650 IATA/ICAO.

Both certificates issued by third-party authorized laboratories.

Applications

For those clients who understand that a rigid secondary container offers added value.

This TRIPLE PACKAGING system for Category B biological samples is recommended for those services in which disposable containers are used.

Advantages

- There are no extra costs for the returning or sterilization of the container, or storage costs in the receiving laboratory.
- 2. It is fully recyclable as clean plastic (unless it has been contaminated).
- All our packaging systems can be shipped via any road or air transport company without needing to comply with any other requirements.

MEDIMAIL TRIPLE PACKAGING:










Primary receptacles are not included.

MEDIMAIL CANISTER MCM1 SALE UNIT:





Carton 30 units. Weight +/-11,5 Kg.

Weight of CANISTER MCM1: 340 g.

Height: 88 cm. Length: 55 cm. Width: 40 cm.

36208 Vigo - Pontevedra (Spain)