



Unique range of In vitro Positive Controls

- ✓ Verification of immunoassay devices
- ✓ Quality management system
- ✓ Education purposes
- ✓ Laboratory personnel training

OXA-48

KPC

NDM


VIM

IMP



P-10R11

RESIST penta O.K.N.V.I control



Innovative solutions for effective diagnostics

Discover the whole Positive Controls range

Reference	Product name	Antigen	To use with Coris BioConcept's kits	Packaging
P-10R11	RESIST penta O.K.N.V.I control (5 vials of individual controls)	OXA-48 KPC NDM VIM IMP	OXA-48 K-SeT (K-15R1) IMP K-SeT (K-15R10) RESIST-3 O.O.K. K-SeT (K-15R4) RESIST-3 O.K.N. K-SeT (K-15R5) O.K.N.V.I. RESIST-5 (K-15R11)* *: use only LY-A buffer to dilute the positive control	Freeze-dried
P-10R4-1	OXA-163 Positive Control	OXA-163	RESIST-3 O.O.K. K-SeT (K-15R4)	Freeze-dried
P-10R7	OXA-23 Positive Control	OXA-23	OXA-23 K-SeT (K-15R7)	Freeze-dried